


ペプチドおよびタンパク質 (MW 4,300 - 25,700)

Peptides and proteins (MW 4,300 - 25,700)

H140219D

- | | |
|-----------------------------------|-------------|
| 1. Cytochrome c | (MW 12,400) |
| 2. Insulin | (MW 5,700) |
| 3. Amyloid β protein (1-40) | (MW 4,300) |
| 4. Lysozyme | (MW 14,000) |
| 5. α -Chymotrypsinogen A | (MW 25,700) |



Column	: Meteoric Core C18 BIO (2.7 μ m, 16 nm) 150 X 2.1 mm I.D.
Eluent	: A) water/TFA (100/0.1) B) acetonitrile/TFA (100/0.1) 25-70%B (0-15 min), 70%B (15-17 min)
Flow rate	: 0.2 mL/min
Temperature	: 40°C
Detection	: UV at 220 nm
Injection	: 2 μ L (0.05-0.2 mg/mL)
Pressure	: 12.8-16.1 MPa (1860-2330 psi)